



# Boston Drug Analysis Laboratory

Massachusetts Department of Public Health  
William Hinton State Laboratory Institute  
305 South Street,  
Jamaica Plain, MA. 02130  
617-983-6622 (Office)  
617-983-6625 (Fax)

# Welcome to the Boston Drug Lab!

## ◆ Mission Statement:

- Identifies unknown substances for local, state and federal agencies in Massachusetts.
- Provides accurate and timely analysis.
- Provide Certificates of Analysis.




# Melendez-Diaz

◆ Trials scheduled: 3065

◆ # of chemist's appearances: 447

Impact of Melendez-Diaz from June 26<sup>th</sup>, 2009 through  
March 31<sup>st</sup>, 2011 on the Boston Drug Analysis Laboratory.

# Evidence Submission



Annie Dookhan

- ◆ At the evidence window, the evidence officer will receive various types of evidence from police departments.
- ◆ Consists of many different types:
  - Substance (Cocaine, Heroin, etc)
  - Rx
  - Marijuana

# Evidence Submission

- ◆ A Drug Receipt is filled out. It will contain the following information:

- Submitting PD and officer.
- Defendant's name (if known)
- A description of the item(s).
- An initial gross weight of the evidence .
- Then each piece of evidence is assigned a laboratory number.

PLEASE PRINT CLEARLY OR TYPE ALL INFORMATION

The Commonwealth of Massachusetts  
Executive Office of Public Safety  
Massachusetts State Police Laboratory Division

**DRUG RECEIPT**

Massachusetts State Police Laboratory  
100 State Street  
Boston, MA 02109  
Phone: 617-725-1000  
Fax: 617-725-1001

City or Department: \_\_\_\_\_ Police Reference No.: \_\_\_\_\_

Name and Rank of Submitting Officer: \_\_\_\_\_  
(Print Name, Rank, Initials)

To be completed by Submitting  
Description of Item Submitted

To be completed by Lab Personnel  
Gross Weight Lab Number


Received by: \_\_\_\_\_ Date: \_\_\_\_\_

# Evidence Submission

The image shows a white evidence submission form with handwritten information. At the top right, a bar code is labeled 'B11-02486'. The form is divided into two main sections. The left section contains fields for 'COHASSET PD', 'Officer' (DEA GORDON HUNN), 'Arrest Date' (3/4/11), 'Agency' (025010101), 'Court', 'Judge' (11-34 35 36 AC), and 'Case #'. The right section contains fields for 'Lab #', 'Date Rec'd', 'POW' (37-11), 'RES', 'VM' (checked), 'Rx', and 'Other'. At the bottom center, there is a circular seal with the word 'EVIDENCE' above it.

- ◆ The sample is received in a heat sealed and initialed plastic bag.
- ◆ The sample is placed into the manila envelope.
- ◆ A bar code with the appropriate Lab # will be adhered to a manila envelope.

# Evidence Submission



No.                      Date Analyzed:

City:

Defendant:

Amt:        Subst:                      Cont:

Date Rec:                      From:

Gross Wt:                      Number of Tests:

Net Weight:

- ◆ An Evidence Control Card is generated based on the information on the drug receipt.
- ◆ The control card is placed in the manila envelope and put into the safe.
- ◆ At this point , the chain of custody is with the Evidence Office



# The Safe



Evidence is transferred to the safe for storage



Annie Dookhan




# Case Assignments



- ◆ An Evidence Officer assigns the samples to individual chemists.
- ◆ The chemist will review the envelope, control card and the sample for accuracy. The chemist accepts custody of the sample by signing the corresponding Log Book.
- ◆ Now, the chain of custody is with the Chemist

# Analysis



- ◆ The Chemist will analyze the substances according to the Lab Policies and Procedures.
- ◆ The Lab follows the recommendation set forth by the Scientific Working Group for the Analysis of Seized Drugs (SWGDRUG).
- ◆ There are different procedures based on the type of sample being analyzed (i.e. a suspected cocaine sample vs. a suspected marijuana sample)

- [illegible]

१०५५

# Preliminary Analysis - Substance

## Weight

- ◆ The weight of the sample can be determined by two methods:


- Indirect
- Direct

- ◆ Sampling Plan

- 10%
- Square Root
- Hypergeometric



# Preliminary Analysis – Heroin




Annie Dookhan

# Preliminary Analysis – Heroin

## Spot Test

- ◆ Spot tests are performed on the substance.
  - ✓ A total of four spot tests
    - Cobalt thiocyanate
    - Marquis
    - Froehdes
    - Meckes
  - ✓ A purple color change for the Marquis and Froehdes, and a green color change for the Meckes, indicates that the substance may be presumptively positive for an opiate.





# Preliminary Analysis – Cocaine




Annie Dookhan

# Preliminary Analysis – Cocaine

## Spot Test

Cobalt Thiocyanate    Marquis    Froehdes    Meckes





◆ Spot tests are performed on the substance.

- ✓ A total of four spot tests
  - Cobalt thiocyanate
  - Marquis
  - Froehdes
  - Meckes
- ✓ A blue color change in the cobalt thiocyanate indicates that substance may be presumptively positive for cocaine.

# Preliminary Analysis – Cocaine

## Microcrystalline Tests




Annie Dookhan

- ◆ Performed if there is a positive cobalt thiocyanate test .
- ◆ Two microcrystalline tests (Structural ID)
  - Gold chloride (AuCl)
  - di-p-toluoyl-L-tartaric acid (TLTA)
- ◆ When these reagents are added to cocaine samples , they will form specific types crystals that can be viewed under a microscope.

# Confirmatory Analysis

Gas Chromatography/Mass Spectrometry (GC-MS)




Annie Dookhan

18

# Confirmatory Analysis

## Gas Chromatography/Mass Spectrometry (GC/MS)

- ◆ For substances from Class A to D, samples are submitted to GC/MS for confirmatory testing.
- ◆ A small amount of the substance is placed in a vial and dissolved with a solvent (methanol).
- ◆ An aliquot of the sample is used when injected into the GC/MS.
- ◆ The instrument separates the components of the substance and then breaks it down into ions, and a chromatogram is generated.
- ◆ The Confirmatory Chemist analyzes the chromatogram of the unknown sample against a known standard to determine the presence of the illicit narcotic.





# Analysis – Marijuana




Annie Dookhan



# Analysis – Marijuana

## Macroscopic identification





- ◆ A macroscopic examination of the sample is done.
- ◆ Identify serrated leaves and morphological characteristics.
- ◆ Does it look like green, leafy, vegetable matter?

# Analysis – Marijuana

## Microscopic identification


- ◆ A microscopic identification of the sample is performed using a stereoscope.
- ◆ Identify two particular types of hairs. (Structural ID)
  - Cystolithic (bear claws)
  - Glandular



Cystolithic

# Analysis – Marijuana

## Modified Duquenois Levine Test



- ◆ A Modified Duquenois Levine test is performed.
  - QC of Reagents
  - The chemist will perform a series of extractions during this test
  - The transfer of purple color from one layer to another, indicates the presence of THC


# Analysis – Rx



Annie Dookhan

# Analysis – Rx

- ◆ Examine & document each pill for identifying marks (shape, color, imprint).
- ◆ Number of pills and the net weight is documented.
- ◆ Compare the ballistics to those in a pharmaceutical reference .
  - Micromedex Software (licensed)
  - Drug Identification Bible
  - Physician's Desk Reference (PDR)
- ◆ Those substances containing an active ingredient from Class A to D will undergo confirmatory testing by GC/MS.
- ◆ Class E substances are confirmed by appearance and labeling .





# Evidence Control Card

- ◆ The evidence control card is filled out by the chemist/s.
  - Results
  - Net Weight
  - Date analyzed
  - Chemist/s who performed the analysis
- ◆ Card is returned to the evidence office.
- ◆ A Certificate of Analysis is generated based on the information on the control card.

No. B11-01559	Date Analyzed: 4-13-11
City: Quincy Police Dept.	
Officer: Detective BRIAN COEN	
Def: MUTTS, JAMES	
Amount:	Subst: SUB
No. Cont: 1	Cont: paper fold
Date Rec'd: 02/09/2011	No. Analyzed:
Gross Wt.: 2.93	Net Weight: 0.05
	# Tests: 6AD
	DDJR
Prelim: Cocaine ©	Findings: Cocaine



# Certificate of Analysis

**The Commonwealth of Massachusetts**  
Executive Office of Health and Human Services  
Department of Public Health  
William A. Hinton State Laboratory Institute  
305 South Street, Jamaica Plain, MA 02130

CEVALA, PATRICIA  
DEPUTY SECRETARY  
UNIVERSITY MEDICAL CENTER  
JENNIFER BUCKLEY, MD  
MEDICAL STAFF  
JENNIFER BUCKLEY, MD  
MEDICAL STAFF

NO. E10-09477  
I hereby certify that the 4 Tablets  
Contained in 1 Plastic Bag  
Submitted by Detective KATHLEEN MULLING of the Needham Police Dept.  
MARKED: E10-09477

Has been examined with the following results:  
The tablets were found to contain:  
Carbamazepine, as defined in Chapter 94C, Controlled Substance Act, Section 1A,  
Class C.

DEPENDANT: PIERCE, STEVEN

ASSISTANT ANALYST  
Annie Dookhan Kate Corbett

On this 07<sup>th</sup> day of December 2010, before me, the undersigned notary public, personally appeared the above signed  
subscriber(s), having proved to me through Department of Public Health documentation to be the person(s) whose name(s)  
appear on this certificate and to be (an) assistant analyst(s) of the Department of Public Health, and who swore to me that  
the contents of this document are truthful and accurate to the best of his/her knowledge and belief.

Daniela Frasca  
Notary Public  
Commonwealth of Massachusetts  
My commission expires on:  
December 17, 2015

Daniela Frasca, NOTARY PUBLIC  
My commission expires on December 17, 2015

Chapter 121, Section 13 of the General Laws  
This certificate shall be sworn to before a Justice of the Peace or Notary Public, and the just shall contain a statement that the  
analyst or assistant analyst of the department, when properly executed, it shall be prima facie evidence of its  
composition, quality, and the net weight of the substance or tablet, drug, poison, medicine, or chemical analyzed, and the court  
shall take notice of the signature of the analyst or assistant analyst, and of the fact that he/she is such.

- ◆ A Certificate of Analysis is generated for each Lab#. It is reviewed and signed by the chemist/s.

- ◆ The first signature is the custodial chemist, who maintains custody of the sample and performs the preliminary analysis & net weight. The second signature is the chemist who performs confirmatory testing.

- ◆ Then the Certificate is notarized.

Annie Dookhan

# Evidence Returns



- ◆ Upon completion, the sample and the Certificate of Analysis is placed in the manila envelop. Both the envelope and the control card will be returned to the safe.
- ◆ Evidence and drug certificates are returned to the appropriate law enforcement agency.

# Court

## Discovery Packets

- ◆ ADAs may request, in writing, a discovery packet from the lab via email or fax.
- ◆ Include the Lab # on the discovery request.
- ◆ 14 days advance notice for discovery packets.
- ◆ Discovery packets will consist of:
  - Drug Receipt
  - Evidence Control Card
  - Chemist Drug Analysis Sheets
  - GC/MS Analytical Data

# Court

## Summons

- ◆ Summons can be emailed, faxed or mailed to the Lab.
- ◆ Include the Lab # on the summons.
- ◆ 30 days advance notice for summons.
- ◆ 48 hours verbal or email follow-up, if needed to appear for court.